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Institute for Molecular Science, Instrument Center

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Liquid Helium (LHe) is a cryogen with a boiling point of 4.2 K (-269 °C).
As handling Liquid Helium is more difficult than Liquid Nitrogen,
users must become proficient in its operation.

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Improper use can lead to serious accidents, including fatal injuries and
equipment damage, affecting not only the user but also the surroundings.
Even if an accident does not occur, failed experiments result
in a significant waste of valuable machine time.

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
Helium is a rare and expensive gas, supplied on the premise that it will be
recycled. To provide the best possible experimental environment within
a limited budget, we ask for your cooperation in improving the recovery rate.

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Please work safely.

If you have any concerns regarding the operation, please consult the staff.

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When liquid vaporizes, it expands 700 times. Gas helium expands significantly as it warms, resulting in burst of sealed containers.

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Secure the gas flow path by pressure relief valves of the dewars and cryostats.

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Touching objects cooled by liquid helium can cause severe frostbite, resulting in blistering.

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Nitrogen and Oxygen in the air freeze at the extremely low temperature. It can block the vessels and cryostats and burst at high pressure, which occurs in much shorter time than by liquid nitrogen exposure.

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Ventilate adequately in a room where a large amount of cryogen vaporizes. Monitor the oxygen concentration ratio.

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Air remaining in the annular space is trapped by the low temperature of liquid helium, creating the high vacuum.

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
**Leave at least 5L of liquid helium to
maintain the high vacuum.**

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Thin layered insulations keep the low temperature, preventing from external thermal radiation.

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The inner tank is suspended by a fragile neck tube made of epoxy resin.

Avoid any shocks or vibrations.

容器のバルブの詳細 / Description of Valves on LHe Vessel (Dewar)

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It releases gas to the atmosphere when the internal pressure exceeds 0.5 bar, then closes automatically.

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A transfer tube is inserted here. Keep it sealed with a plug when not in use.

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Used for filling at the Instrument Center. Keep this valve closed during normal operation.

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Valve for connecting to the recovery line. Compatible with standard recovery hoses used in IMS.

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
It is mainly used for connecting to a balloon or pure helium gas cylinder for pressurization.

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Check the pressure frequently and Keep under 0.5 bar throughout the transfer.

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Gas is vented when the pressure exceeds 1.0-1.5 bar.

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Improper handling of the valves results in air flowing into the vessels and possible accidents. Always verify that valves are properly opened or closed when handling hoses.

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It is a typical vacuum-tight seal for cryostats, dewars, superconducting magnets. If the O-ring is damaged, replace it with a new one.

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
A taper-cut ring compresses O-ring to prevent leaks.

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Do not move the stick without loosening the cover; the O-ring may get dragged into the gap and break the seal.

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
It prevents leakage even when the transfer tube is moving.

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The purity of helium gas in recovery piping is monitored. If the purity is low, locate and seal the leak immediately.

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
The alarm beeps when the oxygen concentration is low. Ventilate or open the door when you use cryogen.

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8L of helium are lost every 100L of use and must be purchased from gas suppliers for replenishment.

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Helium will be more expensive and difficult to get in the future. Please cooperate in the recovery.

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
URL for Helium Recovery Purity Monitor

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各研究室
Laboratories

ヘリウム回収配管
Recovery Piping
(Low-pressure than atmosphere)

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The recovery pipes are -0.2 kPa, lower than atmospheric pressure.
Do not leave the pipes open to the atmosphere.

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
Do not close any valves !

- ②ヘリウム容器を重量計に載せる
Put the liquid helium vessel on the weight scale.

- ③持出返却システムで「持出」を押す
Tap "TakeOut" on the Takeout/Return System.

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




Purge the transfer tube with helium gas before it cools down.

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To prevent the transfer tube from being blocked by frozen air, open the flow control valve of the transfer tube before it is inserted into a liquid helium vessel. Then the air is driven out by the gas helium flow.

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Gaseous helium is a good coolant because of its large sensible heat.

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If 100g of copper pipe at room temperature is cooled by liquid helium → 3 liters vapor

if the pipe is cooled by gas helium in vessel → only 0.08 liters vapor

You can save evaporation of liquid helium when you slowly lower the transfer tube into a cold vessel or cryostat.

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Raise the transfer tube by 10cm up from the bottom of the vessel.

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It prevents the transfer tube from plugging by impurity snow piled on the bottom of the vessel, and makes sure several liters of liquid helium remain in a vessel after the operation.

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Methods of liquid helium transfer depend on the types of cryostats and transfer tubes.

Please consult the staff in charge for details.

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Immediately warm the affected area with running water (below 40°C), then treat it as a normal burn.
Seek medical attention for serious injuries.

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Oxygen levels are low. The alarm is sounding.

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Do not approach or enter the area. Ventilate until oxygen levels return to normal.

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A liquid helium vessel is blocked.

(e.g., cannot insert transfer tube or pressurize, etc.)

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Melt the ice at the block with a warm copper pipe. Never look directly into the vessel port.
Immediately bring the vessel to the cryogenic room A105, or Call the cryogenic staff.

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In case of emergency, refer to the "Emergency Contact Guide" located near the extension phones in each room.

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The liquid helium vessel is completely empty.

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Return the vessel to the cryogenic room A105 as soon as possible.

We will inspect and prepare it for reuse after checking the internal condition.

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The side (not top) of the vessel is cold and covered with water drops.

Evaporation of liquid helium is larger than usual.

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The vacuum insulation may have deteriorated. Please inform the cryogenic staff when you return it.

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The vessel was left open to the atmosphere and may be blocked.

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If in following cases, please contact the cryogenic staff immediately, as the vessel is blocked :

- The vessel pressure does not drop even when the vent valves are open.
- The transfer tube cannot be inserted or lowered to the bottom of the vessel.
- Helium gas for pressurizing cannot be introduced into the vessel.

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Liquid helium cannot be transferred to the cryostat or tank at all.

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The transfer tube or cryostat might be blocked by ice, or Liquid helium vessel might be empty.

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Remove the transfer tube from the cryostat and vessel, then warm and dry the inside with dry gas.

To prevent plugging, purge the air in the transfer tube with helium gas before inserting into the vessel.

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The insulated part of the transfer tube is cold, and covered in frost or water drops.

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The vacuum has deteriorated and needs to be re-evacuated. Please consult the equipment staff.

You can use it as it is, but the efficiency of transfer should be bad.

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Pressure of the vessel does not decrease even if the vent valve is open when transferring.

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Verify the recovery valve on the lab's wall is open, and the hoses are not kinked or crushed.

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When the purity of the helium gas recovery line decreases, please locate and seal the leak immediately.
If locating the leak is difficult, test with a handy leak detector which is available in the cryogenic room.

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Myodaiji Cryogenic Staff, Instrument Center

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Yamate Cryogenic Staff, Instrument Center

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